



2 X  $\phi$ .119 [3]  $\nabla$ .16 [4]

3 X 1/4-20 UNC  $\nabla$ .3 [7.6]

2 X M5X.8  $\nabla$ .3 [7.6]

UNLESS OTHERWISE SPECIFIED:		NAME	DATE	<b>RESONON, INC.</b>	
DIMENSIONS ARE IN INCHES TOLERANCES: FRACTIONAL $\pm 1/64$ ANGULAR: MACH $\pm 0.1$ BEND $\pm 1$ TWO PLACE DECIMAL $\pm 0.01$ THREE PLACE DECIMAL $\pm 0.005$ BREAK ALL SHARP EDGES		DRAWN: Dodge	06/19/2019	123 Commercial Drive, Bozeman, MT 59715 406.586.3356	
DO NOT SCALE DRAWING		PROJECT:		TITLE:	
MATERIAL N.A.		COMMENTS:		Pika IR-L	
FINISH ROUGHNESS TO BE 32 RMS FOR ALL NON-THREADED SURFACES N.A.		<b>PROPRIETARY AND CONFIDENTIAL</b> THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF RESONON, INC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF RESONON, INC. IS PROHIBITED.		SIZE <b>A</b>	PART NO.
				REV:	SHEET 1 OF 1

Note: Shown without objective lens.